

## Box size 32

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wx dx h mm)
<b>MUF 40</b>	32	2 x USR34-28	2 x 8	32	139 x 278 x 292
<b>ULUF P10</b>	32	-	-	9	-
<b>ULUF P50</b>	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
<b>ULUF P90</b>	32	4 x USR34-38	3 x 8	96	139 x 414 x 292
<b>ULUF P390</b>	32	9 x USR34-48	4 x 8	288	139 x 549 x 292
		3 x USR34-49	4 x 9	108	139 x 549 x 328
		Total: 12 racks	-	(= 396)	-
<b>ULUF P390 GG</b>	32	12 x USR34-48	4 x 8	384	139 x 549 x 256
<b>ULUF P500</b>	32	12 x USR34-48	4 x 8	384	139 x 549 x 292
		4 x USR34-49	4 x 9	144	139 x 549 x 328
		Total: 16 racks	-	(=528)	-
<b>ULUF P500 GG</b>	32	16 x USR34-48	4 x 8	512	139 x 549 x 292
<b>ULUF P610</b>	32	15 x USR34-48	4 x 9	480	139 x 549 x 292
		5 x USR34-49	4 x 9	180	139 x 549 x 328
		Total: 20 racks	-	(=660)	-
<b>ULUF P610 GG</b>	32	20 x USR34-48	4 x 8	640	139 x 549 x 292
<b>ULUF P820</b>	32	21 x USR34-48	4 x 8	672	139 x 549 x 292
		7 x USR34-49	4 x 9	252	139 x 549 x 328
		Total: 28 racks	-	(=924)	-
<b>ULUF P820 GG</b>	32	28 x USR34-48	4 x 8	784	139 x 549 x 292
<b>ULUF 15</b>	32	-	-	9	-
<b>ULUF 65</b>	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
<b>ULUF 60</b>	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
<b>ULUF 120</b>	32	4 x USR34-37	3 x 7	84	139 x 414 x 256
<b>ULUF 125</b>	32	4 x USR34-37	3 x 7	84	139 x 414 x 256
<b>ULUF 400</b>	32	12 x USR34-46	4 x 6	288	139 x 549 x 220
		3 x USR34-48	4 x 8	96	139 x 549 x 292
		Total: 15 racks	-	(=384)	-
<b>ULUF 450-2M</b>	32	15 x USR34-46	4 x 6	360	139 x 549 x 220
<b>ULUF 500</b>	32	16 x USR34-46	4 x 6	384	139 x 549 x 220
		4 x USR34-48	4 x 8	128	139 x 549 x 292
		Total: 20 racks	-	(=512)	-
<b>ULUF 550-2M</b>	32	20 x USR34-46	4 x 6	480	139 x 549 x 220
<b>ULUF 700</b>	32	20 x USR34-46	4 x 6	480	139 x 549 x 220
		5 x USR34-48	4 x 8	160	139 x 549 x 292
		Total: 25 racks	-	(=640)	-
<b>ULUF 750-2M</b>	32	25 x USR34-46	4 x 6	600	139 x 549 x 220

Box size 40

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	40	2 x USR40-26	2 x 6	24	139 x 278 x 256
<b>ULUF P10</b>	40	-	-	7	-
<b>ULUF P50</b>	40	2 x USR40-37	3 x 7	42	139 x 414 x 298
<b>ULUF P90</b>	40	4 x USR40-37	3 x 7	84	139 x 414 x 298
<b>ULUF P390</b>	40	9 x USR40-47	4 x 7	252	139 x 549 x 298
		3 x USR40-48	4 x 8	96	139 x 549 x 340
		Total: 12 racks	-	(=348)	-
<b>ULUF P390 GG</b>	40	12 x USR40-47	4 x 7	336	139 x 549 x 298
<b>ULUF P500</b>	40	12 x USR40-47	4 x 7	336	139 x 549 x 298
		4 x USR40-48	4 x 8	128	139 x 549 x 340
		Total: 16 racks	-	(=464)	-
<b>ULUF P500 GG</b>	40	16 x USR40-47	4 x 7	448	139 x 549 x 298
<b>ULUF P610</b>	40	15 x USR40-47	4 x 7	420	139 x 549 x 298
		5 x USR40-48	4 x 8	160	139 x 549 x 340
		Total: 20 racks	-	(=580)	-
<b>ULUF P610 GG</b>	40	20 x USR40-47	4 x 7	560	139 x 549 x 298
<b>ULUF P820</b>	40	21 x USR40-47	4 x 7	588	139 x 549 x 298
		7 x USR40-48	4 x 8	224	139 x 549 x 340
		Total: 28 racks	-	(=812)	-
<b>ULUF P820 GG</b>	40	28 x USR40-47	4 x 7	588	139 x 549 x 298
<b>ULUF 15</b>	40	-	-	7	-
<b>ULUF 65</b>	40	2 x USR40-37	3 x 7	42	139 x 414 x 298
<b>ULUF 60</b>	40	2 x USR40-37	3 x 7	42	139 x 414 x 298
<b>ULUF 120</b>	40	4 x USR40-36	3 x 6	72	139 x 414 x 256
<b>ULUF 125</b>	40	4 x USR40-36	3 x 6	72	139 x 414 x 256
<b>ULUF 400</b>	40	12 x USR40-45	4 x 5	240	139 x 549 x 214
		3 x USR40-47	4 x 7	84	139 x 549 x 298
		Total: 15 racks	-	(=324)	-
<b>ULUF 450-2M</b>	40	15 x USR40-45	4 x 5	300	139 x 549 x 214
<b>ULUF 500</b>	40	16 x USR40-45	4 x 5	320	139 x 549 x 214
		4 x USR40-47	4 x 7	112	139 x 549 x 298
		Total: 20 racks	-	(=432)	-
<b>ULUF 550-2M</b>	40	20 x USR40-45	4 x 5	400	139 x 549 x 214
<b>ULUF 700</b>	40	20 x USR40-45	4 x 5	400	139 x 549 x 214
		5 x USR40-47	4 x 7	140	139 x 549 x 298
		Total: 25 racks	-	(=540)	-
<b>ULUF 750-2M</b>	40	25 x USR40-45	4 x 5	500	139 x 549 x 214

Box size 50

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	50	2 x USR50-25	2 x 5	20	141 x 286 x 272
<b>ULUF P10</b>	50	-	-	5	-
<b>ULUF P50</b>	50	2 x USR50-35	3 x 5	30	140 x 424 x 274
<b>ULUF P90</b>	50	4 x USR50-35	3 x 5	60	140 x 424 x 274
<b>ULUF P390</b>	50	9 x USR50-45	4 x 5	180	140 x 558 x 274
		3 x USR50-46	4 x 6	72	140 x 558 x 328
		Total: 12 racks	-	(=252)	-
<b>ULUF P390 GG</b>	50	6 x USR50-45	4 x 5	120	140 x 558 x 274
		6 x USR50-46	4 x 6	144	140 x 558 x 328
		Total: 12 racks	-	(= 264)	-
<b>ULUF P500</b>	50	12 x USR50-45	4 x 5	240	140 x 448 x 274
		4 x USR50-46	4 x 6	96	140 x 558 x 328
		Total: 16 racks	-	(=336)	-
<b>ULUF P500 GG</b>	50	8 x USR50-45	4 x 5	160	140 x 448 x 274
		8 x USR50-46	4 x 6	192	140 x 558 x 328
		Total: 16 racks	-	(= 352)	-
<b>ULUF P610</b>	50	15 x USR50-45	4 x 5	300	149 x 558 x 274
		5 x USR50-46	4 x 6	120	140 x 558 x 328
		Total: 20 racks	-	(=420)	-
<b>ULUF P610 GG</b>	50	10 x USR50-45	4 x 5	200	149 x 558 x 274
		10 x USR50-46	4 x 6	240	140 x 558 x 328
		Total: 20 racks	-	(= 440)	-
<b>ULUF P820</b>	50	21 x USR50-45	4 x 5	420	140 x 558 x 274
		7 x USR50-46	4 x 6	168	140 x 558 x 328
		Total: 28 racks	-	(=588)	-
<b>ULUF P820 GG</b>	50	14 x USR50-45	4 x 5	280	140 x 558 x 274
		14 x USR50-46	4 x 6	336	140 x 558 x 328
		Total: 28 racks	-	(= 616)	-
<b>ULUF 15</b>	50	-	-	5	-
<b>ULUF 65</b>	50	2 x USR50-35	3 x 5	30	140 x 424 x 274
<b>ULUF 60</b>	50	2 x USR50-35	3 x 5	30	140 x 424 x 274
<b>ULUF 120</b>	50	4 x USR50-34	3 x 4	48	140 x 424 x 219
<b>ULUF 125</b>	50	4 x USR50-34	3 x 4	48	140 x 424 x 219
<b>ULUF 400</b>	50	12 x USR50-44	4 x 4	192	140 x 566 x 219
		3 x USR50-45	4 x 5	60	140 x 558 x 274
		Total: 15 racks	-	(= 252)	-
<b>ULUF 450-2M</b>	50	15 x USR50-44	4 x 4	240	140 x 566 x 219
<b>ULUF 500</b>	50	16 x USR50-44	4 x 4	256	140 x 566 x 219
		4 x USR50-45	4 x 5	80	140 x 558 x 274
		Total: 20 racks	-	(= 336)	-
<b>ULUF 550-2M</b>	50	20 x USR50-44	4 x 4	320	140 x 566 x 219
<b>ULUF 700</b>	50	20 x USR50-44	4 x 4	320	140 x 566 x 219
		5 x USR50-45	4 x 5	100	140 x 558 x 274
		Total: 25 racks	-	(= 420)	-
<b>ULUF 750-2M</b>	50	25 x USR50-44	4 x 4	200	140 x 566 x 219

Box size 75

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	75	2 x USR75-23	2 x 3	12	140 x 286 x 242
<b>ULUF P10</b>	75	-	-	4	-
<b>ULUF P50</b>	75	2 x USR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF P90</b>	75	4 x USR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF P390</b>	75	9 x USR75-43	4 x 3	108	140 x 558 x 240
		3 x USR75-44	4 x 4	48	140 x 558 x 319
		Total: 12 racks	-	(= 156)	-
<b>ULUF P390 GG</b>	75	6 x USR75-43	4 x 3	72	140 x 558 x 240
		6 x USR75-44	4 x 4	96	140 x 558 x 319
		Total: 12 racks	-	(= 168)	-
<b>ULUF P500</b>	75	12 x USR75-43	4 x 3	144	140 x 558 x 240
		4 x USR75-44	4 x 4	64	140 x 558 x 319
		Total: 16 racks	-	(= 208)	-
<b>ULUF P500 GG</b>	75	8 x USR75-43	4 x 3	96	140 x 558 x 240
		8 x USR75-44	4 x 4	128	140 x 558 x 319
		Total: 16 racks	-	(= 224)	-
<b>ULUF P610</b>	75	15 x USR75-43	4 x 3	180	140 x 558 x 240
		5 x USR75-44	4 x 4	80	140 x 558 x 319
		Total: 20 racks	-	(= 260)	-
<b>ULUF P610 GG</b>	75	10 x USR75-43	4 x 3	120	140 x 558 x 240
		10 x USR75-44	4 x 4	160	140 x 558 x 319
		Total: 20 racks	-	(= 280)	-
<b>ULUF P820</b>	75	21 x USR75-43	4 x 3	252	140 x 558 x 240
		7 x USR75-44	4 x 4	112	140 x 558 x 319
		Total: 28 racks	-	(= 364)	-
<b>ULUF P820 GG</b>	75	14 x USR75-43	4 x 3	168	140 x 558 x 240
		14 x USR75-44	4 x 4	224	140 x 558 x 319
		Total: 28 racks	-	(= 392)	-
<b>ULUF 15</b>	75	-	-	4	-
<b>ULUF 65</b>	75	2 x USR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF 60</b>	75	2 x USR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF 120</b>	75	4 x USR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF 125</b>	75	4 x USR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF 400</b>	75	12 x USR75-43	4 x 3	144	140 x 558 x 240
		3 x USR75-44	4 x 4	192	140 x 558 x 319
		Total: 15 racks	-	(= 192)	-
<b>ULUF 450-2M</b>	75	15 x USR75-43	4 x 3	180	140 x 558 x 240
<b>ULUF 500</b>	75	16 x USR75-43	4 x 3	192	140 x 558 x 240
		4 x USR75-44	4 x 4	64	140 x 558 x 319
		Total: 20 racks	-	(= 256)	-
<b>ULUF 550-2M</b>	75	20 x USR75-43	4 x 3	240	140 x 558 x 240
<b>ULUF 700</b>	75	20 x USR75-43	4 x 3	240	140 x 558 x 240
		5 x USR75-44	4 x 4	80	140 x 558 x 319
		Total: 25 racks	-	(= 320)	-
<b>ULUF 750-2M</b>	75	25 x USR75-43	4 x 3	300	140 x 558 x 240

Box size 80

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	80	2 x USR80-23	2 x 3	12	139 x 278 x 250
<b>ULUF P10</b>	80	-	-	3	-
<b>ULUF P50</b>	80	2 x USR80-33	3 x 3	18	139 x 414 x 250
<b>ULUF P90</b>	80	4 x USR80-33	3 x 3	36	139 x 414 x 250
<b>ULUF P390</b>	80	9 x USR80-43	4 x 3	108	139 x 549 x 250
		3 x USR80-44	4 x 4	48	139 x 549 x 332
		Total: 12 racks	-	(= 156)	-
<b>ULUF P390 GG</b>	80	6 x USR80-43	4 x 3	72	139 x 549 x 250
		6 x USR80-44	4 x 4	96	139 x 549 x 332
		Total: 12 racks	-	(= 168)	-
<b>ULUF P500</b>	80	12 x USR80-43	4 x 3	144	139 x 549 x 250
		4 x USR80-44	4 x 4	64	139 x 549 x 332
		Total: 16 racks	-	(= 208)	-
<b>ULUF P500 GG</b>	80	8 x USR80-43	4 x 3	96	139 x 549 x 250
		8 x USR80-44	4 x 4	128	139 x 549 x 332
		Total: 16 racks	-	(= 224)	-
<b>ULUF P610</b>	80	15 x USR80-43	4 x 3	180	139 x 549 x 250
		5 x USR80-44	4 x 4	80	139 x 549 x 332
		Total: 20 racks	-	(= 260)	-
<b>ULUF P610 GG</b>	80	10 x USR80-43	4 x 3	120	139 x 549 x 250
		10 x USR80-44	4 x 4	160	139 x 549 x 332
		Total: 20 racks	-	(= 280)	-
<b>ULUF P820</b>	80	21 x USR80-43	4 x 3	252	139 x 549 x 250
		7 x USR80-44	4 x 4	112	139 x 549 x 332
		Total: 28 racks	-	(= 364)	-
<b>ULUF P820 GG</b>	80	14 x USR80-43	4 x 3	168	139 x 549 x 250
		14 x USR80-44	4 x 4	224	139 x 549 x 332
		Total: 28 racks	-	(= 392)	-
<b>ULUF 15</b>	80	-	-	3	-
<b>ULUF 65</b>	80	2 x USR80-33	3 x 3	18	139 x 414 x 250
<b>ULUF 60</b>	80	2 x USR80-33	3 x 3	18	139 x 414 x 250
<b>ULUF 120</b>	80	4 x USR80-33	3 x 3	36	139 x 414 x 250
<b>ULUF 125</b>	80	4 x USR80-33	3 x 3	36	139 x 414 x 250
<b>ULUF 400</b>	80	15 x USR80-43	4 x 3	180	139 x 549 x 250
<b>ULUF 450-2M</b>	80	15 x USR80-43	4 x 3	180	139 x 549 x 250
<b>ULUF 500</b>	80	20 x USR80-43	4 x 3	240	139 x 549 x 250
<b>ULUF 550-2M</b>	80	20 x USR80-43	4 x 3	240	139 x 549 x 250
<b>ULUF 700</b>	80	25 x USR80-43	4 x 3	300	139 x 549 x 250
<b>ULUF 750-2M</b>	80	25 x USR80-42	4 x 3	200	139 x 549 x 250

## Box size 90/96/100

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	90/96/100	2 x USR100-22	2 x 2	8	140 x 286 x 211
<b>ULUF P10</b>	90/96/100	-	-	3	-
<b>ULUF P50</b>	90/96/100	2 x USR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF P90</b>	90/96/100	4 x USR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF P390</b>	90/96/100	9 x USR100-42	4 x 2	72	140 x 566 x 212
		3 x USR100-43	4 x 3	36	140 x 566 x 318
		Total: 12 racks	-	(= 108)	-
<b>ULUF P390 GG</b>	90/96/100	6 x USR100-42	4 x 2	48	140 x 566 x 212
		6 x USR100-43	4 x 3	72	140 x 566 x 318
		Total: 12 racks	-	(= 120)	-
<b>ULUF P500</b>	90/96/100	12 x USR100-42	4 x 2	96	140 x 566 x 212
		4 x USR100-43	4 x 3	48	140 x 566 x 318
		Total: 16 racks	-	(= 144)	-
<b>ULUF P500 GG</b>	90/96/100	8 x USR100-42	4 x 2	64	140 x 566 x 212
		8 x USR100-43	4 x 3	96	140 x 566 x 318
		Total: 16 racks	-	(= 160)	-
<b>ULUF P610</b>	90/96/100	15 x USR100-42	4 x 2	120	140 x 566 x 212
		5 x USR100-43	4 x 3	60	140 x 566 x 318
		Total: 20 racks	-	(= 180)	-
<b>ULUF P610 GG</b>	90/96/100	10 x USR100-42	4 x 2	80	140 x 566 x 212
		10 x USR100-43	4 x 3	120	140 x 566 x 318
		Total: 20 racks	-	(= 200)	-
<b>ULUF P820</b>	90/96/100	21 x USR100-42	4 x 2	168	140 x 566 x 212
		7 x USR100-43	4 x 3	84	140 x 566 x 318
		Total: 28 racks	-	(= 252)	-
<b>ULUF P820 GG</b>	90/96/100	14 x USR100-42	4 x 2	112	140 x 566 x 212
		14 x USR100-43	4 x 3	168	140 x 566 x 318
		Total: 28 racks	-	(= 280)	-
<b>ULUF 15</b>	90/96/100	-	-	3	-
<b>ULUF 65</b>	90/96/100	2 x USR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF 60</b>	90/96/100	2 x USR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF 120</b>	90/96/100	4 x USR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF 125</b>	90/96/100	4 x USR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF 400</b>	90/96/100	12 x USR100-42	4 x 2	96	140 x 566 x 212
		3 x USR100-43	4 x 3	36	140 x 566 x 318
		Total: 15 racks	-	(= 132)	-
<b>ULUF 450-2M</b>	90/96/100	15 x USR100-42	4 x 2	120	140 x 566 x 212
<b>ULUF 500</b>	90/96/100	16 x USR100-42	4 x 2	128	140 x 566 x 212
		4 x USR100-43	4 x 3	48	140 x 566 x 318
		Total: 20 racks	-	(= 176)	-
<b>ULUF 550-2M</b>	90/96/100	20 x USR100-42	4 x 2	160	140 x 566 x 212
<b>ULUF 700</b>	90/96/100	20 x USR100-42	4 x 2	160	140 x 566 x 212
		5 x USR100-43	4 x 3	60	140 x 566 x 318
		Total: 25 racks	-	(= 220)	-
<b>ULUF 750-2M</b>	90/96/100	25 x USR100-42	4 x 2	200	140 x 566 x 212