



#	Description	Comment
1	Protective earth connection	
2	Electrical box	
3	Condensing unit base	
4	Compressor	NLV 12,6CN
5	Inverter	105N4760
6		
7		
8	Compressor sensor	NTC
9		
10	Chamber sensor	PT100
11	Monitoring sensor	PT100
12	Door switch	MAM-1812 B-1
13	Battery	GPL1272 F2FR
14	Power supply	110-230VAC/12VDC LRS-75-12
15	Ferrite	
16	Canbus wire	
17		
18	Condenser fan	NIQ3212 ebmPapst
19	Frame heater	
20	Remote alarm	
A	3 x 1,5mm ²	
B	4 x 1,0mm ²	(1-L, 2-N)
C	3 x 1,5mm ²	
D	4 x 1,0mm ²	(1-L, 2-N)
E	2 x 1,0mm ²	
F	1 x 2,0mm ²	
G	1 x 2,0mm ²	
H	4 x 1,0mm ²	(1-L, 2-N, 3-comp)
I	2 x 1,0mm ²	

Description of changes:		Initials:	Date:
Author: <i>S. S. Jacobs</i>	Checked:	Title: Wiring diagram for ULUF P65MV	
Date: 2024-08-12	Date:	Drawing number: 5086201	
Format: A3	Weight: n/a	Rev: 01	Rev. Date: 12/8-24
Scale: n/a	Weight: n/a	Sheet: 1/1	

All rights to this drawing belong to Arctiko A/S. This drawing may not be duplicated without written permission from Arctiko A/S