



#	Description	Comment
1	Protective earth connection	
2	Electrical box	
3	Condensing unit base	
4	Compressor #1	NLV 8CN
5	Inverter #1	105N4760/105N4761
6	Compressor #2	NLV 8CN
7	Inverter #2	105N4760/105N4761
8	Compressor #1 sensor	NTC
9	Compressor #2 sensor	NTC
10	Chamber sensor	PT100
11	Monitoring sensor	PT100
12	Door switch	MAM-1812 B-1
13	Battery	GPL1272 F2FR
14	Power supply	110-230VAC/12VDC
		LRS-75-12
15	Ferrite	
16	Canbus wire	
17	Lightning	LED
18	Condenser fan	
19	Vacuum valve heater	NIQ3212 ebmPapst
20	Remote alarm	
A	3 x 1,5mm ²	
B		
C	3 x 1,5mm ²	
D	4 x 1,0mm ²	(1-L, 2-N)
E	2 x 1,0mm ²	
F	1 x 2,0mm ²	
G	1 x 2,0mm ²	
H	4 x 1,0mm ²	(1-L, 2-N, 3-comp)
I		
J	2 x 1,0mm ²	

POWER SUPPLY
115- 230 V 50/60 Hz

Description of changes:		Initials:	Date:
Author: <i>S. S. Stubb</i>	Checked:	Title: Wiring diagram for ULUF P400MV	
Date: 2024-08-28	Date:	Drawing number: 5086202	
Format: A3	Weight: n/a	Rev: 01	Rev. Date: 28/8-24
Scale: n/a	Weight: n/a	Sheet: 1/1	

ARCTIKO
true cooling specialists
Lillebæltsvej 90
6715 Esbjerg N. Denmark
www.arctiko.com

All rights to this drawing belong to Arctiko A/S. This drawing may not be duplicated without written permission from Arctiko A/S