


4 pcs plastic coated shelves
max. load kg.30



remote alarm

Description of changes:		Initials:	Date:
Author: <i>S. Stucki</i>	Checked:	Title: PR 350	
Date: 2025-04-28	Date:	Drawing number: 5088265	
Format: A3		Drawing number: 5088265	
Scale: 1:15	Weight: 110kg	Rev: 01	Rev. Date: 28/4-25 Sheet: 1/1
 true cooling specialists Lillebæltsvej 90 6715 Esbjerg N. Denmark www.arctiko.com			