

# PERFORMA-THERM™ MEDICAL TRANSPORT COOLER

One Cooler. Three Configurations. Zero Compromise.

## A next-generation medical transport cooler designed for unmatched performance.

From refrigerated 2–8 °C for vaccines and biologics, to controlled room temperature at 18 °C for pharmaceuticals and diagnostic kits, and frozen storage at –20 °C for plasma and specialty medicines.

This next-generation medical transport cooler gives you the flexibility to handle it all. Built for reliability and confidence in every journey, it ensures your critical payloads arrive exactly as intended, without compromise.


### Applications

Blood, plasma, &  
organ transport

Vaccines & biologics

Clinical trial materials

Hospital, lab, &  
field operations



PERFORMA-THERM™  
MEDICAL TRANSPORT COOLER



### Three Configurations

- Includes 2–8 °C PCM inserts
- Optional -20 °C and Controlled Room Temperature (CRT) panels available
- Swap panels to match your shipment—no dry ice, no mess

### Smarter Portability

- Ergonomic straps and multiple carry options
- Exterior pockets for data loggers, manifests, and documents
- Lightweight < 10 lb.

### Easy to Clean

- Removable hard plastic liner
- Wipe down with alcohol for rapid re-sterilization

### Secure Transport

- Double zipper + double buckle closure system
- Keeps contents safe and compliant

### Stability You Can Count On

- Non-slip rubber feet provide traction.
- Prevents sliding or tipping during transport.

Switch to Smarter Cold Transport with Performa-Therm™  
973-300-9100 Sales Dept. | [thermcoproducts.com/performatherm](http://thermcoproducts.com/performatherm)

## PCM Panels:

State-of-the-art passive technologies  
PCM & VIP extend cold chain  
protection without power



### ACCPTPC

(Replacement PCM panels,  
Set of six (6), 6" x 6" x 1", includes  
qualification, interlocking)



### ACCPTFRZ

(Replacement PCM panels,  
Freezer (-21C), Set of six (6),  
6" x 6" x 1", includes qualification,  
interlocking)



### ACCPTCRT

(Replacement PCM panels,  
Controlled Room Temperature (18C)  
Set of six (6), 6" x 6" x 1", includes  
qualification, interlocking)

## Liner:

Reduce risk of  
contamination  
and leakage



## Outershell:

Heavy duty buckle clips, grip feet,  
and double zipper closure ensure  
proper and safe sealing of the cooler



# ACCPT3L & ACCPT5L Specs



ISTA Certified 7D Duration Results. This product has been tested against both the ISTA 7D & 7E Summer and Winter profiles.

Catalog No.	Payload Volume (L)	Payload Dimensions (in)	Outer Dimensions (in)	Tare Weight (lbs)	Temp Range (c)	Season ISTA 7D	Hours	Temp Range (C)	Season ISTA 7D	Hours	Temp Range (C)	Season ISTA 7E	Hours
ACCPT3L	3	5 x 5 x 5	11.81 x 10.24 x 11.02	9.3	2° to 8° (C)	Summer	48	15° to 25° (C)	Summer	95	-15° to -25° (C)	Summer	43
						Winter	144+		Winter	46		Winter	60
ACCPT5L	5	9 x 5 x 5	15.9 x 10.4 x 10.4	12	2° to 8° (C)	Summer	48	15° to 25° (C)	Summer	96	-15° to -25° (C)	Summer	40
						Winter	144+		Winter	101		Winter	65

## Validated Performance

- ISTA 7E/7D S/W performance tested
- Vacuum Insulated Panel (VIP) Temp-Shield™ system
- Accurate & stable temperature maintaining profile  
2°- 8°C, 18°C or -20 °C for 2+ days

## Routine Vaccine Transport:

[www.cdc.gov/vaccines/hcp/admin/storage/toolkit/storage-handling-toolkit.pdf](http://www.cdc.gov/vaccines/hcp/admin/storage/toolkit/storage-handling-toolkit.pdf)



## Data Loggers:

Ensure the integrity of temperature-sensitive payload with our CIMSense family of digital data loggers. Delivering reliable, continuous cold-chain temperature monitoring precise and excursion tracking throughout transportation.

Trusted names include: LogTag, Lascar and AccuTherm.

*Data Loggers Sold Separately\**

973-300-9100 Sales Dept. | [thermcoproducts.com/performatherm](http://thermcoproducts.com/performatherm)