

CRYO BOX 10x10 / 9x9

SPECIFICATION	PBN05	PBN10	PBN20	PBN40	PBN50	PBN50	PBN75	PBN96
Item no.	DAI 3160	DAI 3161	DAI 3162	DAI 3163	DAI 3164	DAI 3100	DAI 3106	DAI 3112
Box size (height mm)	32	40	50	80	90	50	75	96
Dimensions (WxDxH mm)	130x130x32	130x130x40	130x130x50	130x130x80	130x130x90	133x133x50	133x133x75	133x133x96
Total vial capacity	100 pcs.	100 pcs.	100 pcs.	100 pcs.	100 pcs.	81 pcs.	81 pcs.	81 pcs.
Minimum order quantity	1 pcs. (10x10 box)	1 pcs. (10x10 box)	1 pcs. (10x10 box)	1 pcs. (10x10 box)	1 pcs. (10x10 box)	1 pcs. (9x9 box)	1 pcs. (9x9 box)	1 pcs. (9x9 box)