

Box size 32

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	32	-	-	9	-
ULUF P50MV	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
ULUF P90MV	32	4 x USR34-38	3 x 8	96	139 x 414 x 292
ULUF P400MV	32	12 x USR34-48	4 x 8	384	139 x 549 x 256
ULUF P510MV	32	16 x USR34-48	4 x 8	512	139 x 549 x 292
ULUF P620MV	32	20 x USR34-48	4 x 8	640	139 x 549 x 292
ULUF P830MV	32	28 x USR34-48	4 x 8	784	139 x 549 x 292
ULUF 450-2M	32	15 x USR34-46	4 x 6	360	139 x 549 x 220
ULUF 750-2M	32	25 x USR34-46	4 x 6	600	139 x 549 x 220

Box size 40

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	40	-	-	7	-
ULUF P50MV	40	2 x USR40-37	3 x 7	42	139 x 414 x 298
ULUF P90MV	40	4 x USR40-37	3 x 7	84	139 x 414 x 298
ULUF P400MV	40	12 x USR40-47	4 x 7	336	139 x 549 x 298
ULUF P510MV	40	16 x USR40-47	4 x 7	448	139 x 549 x 298
ULUF P620MV	40	20 x USR40-47	4 x 7	560	139 x 549 x 298
ULUF P830MV	40	28 x USR40-47	4 x 7	588	139 x 549 x 298
ULUF 450-2M	40	15 x USR40-45	4 x 5	300	139 x 549 x 214
ULUF 750-2M	40	25 x USR40-45	4 x 5	500	139 x 549 x 214

Box size 50

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wx dxh mm)
ULUF P10MV	50	-	-	5	-
ULUF P50MV	50	2 x USR50-35	3 x 5	30	140 x 424 x 274
ULUF P90MV	50	4 x USR50-35	3 x 5	60	140 x 424 x 274
ULUF P400MV	50	6 x USR50-45	4 x 5	120	140 x 558 x 274
		6 x USR50-46	4 x 6	144	140 x 558 x 328
		Total: 12 racks	-	(= 264)	-
ULUF P510MV	50	8 x USR50-45	4 x 5	160	140 x 448 x 274
		8 x USR50-46	4 x 6	192	140 x 558 x 328
		Total: 16 racks	-	(= 352)	-
ULUF P620MV	50	10 x USR50-45	4 x 5	200	149 x 558 x 274
		10 x USR50-46	4 x 6	240	140 x 558 x 328
		Total: 20 racks	-	(= 440)	-
ULUF P830MV	50	14 x USR50-45	4 x 5	280	140 x 558 x 274
		14 x USR50-46	4 x 6	336	140 x 558 x 328
		Total: 28 racks	-	(= 616)	-
ULUF 450-2M	50	15 x USR50-44	4 x 4	240	140 x 566 x 219
ULUF 750-2M	50	25 x USR50-44	4 x 4	200	140 x 566 x 219

Box size 75

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	75	-	-	4	-
ULUF P50MV	75	2 x USR75-33	3 x 3	18	140 x 420 x 242
ULUF P90MV	75	4 x USR75-33	3 x 3	36	140 x 420 x 242
ULUF P400MV	75	6 x USR75-43	4 x 3	72	140 x 558 x 240
		6 x USR75-44	4 x 4	96	140 x 558 x 319
		Total: 12 racks	-	(= 168)	-
ULUF P510MV	75	8 x USR75-43	4 x 3	96	140 x 558 240
		8 x USR75-44	4 x 4	128	140 x 558 x 319
		Total: 16 racks	-	(= 224)	-
ULUF P620MV	75	10 x USR75-43	4 x 3	120	140 x 558 x 240
		10 x USR75-44	4 x 4	160	140 x 558 x 319
		Total: 20 racks	-	(= 280)	-
ULUF P830MV	75	14 x USR75-43	4 x 3	168	140 x 558 x 240
		14 x USR75-44	4 x 4	224	140 x 558 x 319
		Total: 28 racks	-	(= 392)	-
ULUF 450-2M	75	15 x USR75-43	4 x 3	180	140 x 558 x 240
ULUF 750-2M	75	25 x USR75-43	4 x 3	300	140 x 558 x 240

Box size 80

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	80	-	-	3	-
ULUF P50MV	80	2 x USR80-33	3 x 3	18	139 x 414 x 250
ULUF P90MV	80	4 x USR80-33	3 x 3	36	139 x 414 x 250
ULUF P400MV	80	6 x USR80-43	4 x 3	72	139 x 549 x 250
		6 x USR80-44	4 x 4	96	139 x 549 x 332
		Total: 12 racks	-	(= 168)	-
ULUF P510MV	80	8 x USR80-43	4 x 3	96	139 x 549 x 250
		8 x USR80-44	4 x 4	128	139 x 549 x 332
		Total: 16 racks	-	(= 224)	-
ULUF P620MV	80	10 x USR80-43	4 x 3	120	139 x 549 x 250
		10 x USR80-44	4 x 4	160	139 x 549 x 332
		Total: 20 racks	-	(= 280)	-
ULUF P830MV	80	14 x USR80-43	4 x 3	168	139 x 549 x 250
		14 x USR80-44	4 x 4	224	139 x 549 x 332
		Total: 28 racks	-	(= 392)	-
ULUF 450-2M	80	15 x USR80-43	4 x 3	180	139 x 549 x 250
ULUF 750-2M	80	25 x USR80-42	4 x 3	200	139 x 549 x 250

Box size 90/96/100

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	90/96/100	-	-	3	-
ULUF P50MV	90/96/100	2 x USR100-32	3 x 2	12	140 x 420 x 212
ULUF P90MV	90/96/100	4 x USR100-32	3 x 2	24	140 x 420 x 212
ULUF P400MV	90/96/100	6 x USR100-42	4 x 2	48	140 x 566 x 212
		6 x USR100-43	4 x 3	72	140 x 566 x 318
		Total: 12 racks	-	(= 120)	-
ULUF P510MV	90/96/100	8 x USR100-42	4 x 2	64	140 x 566 x 212
		8 x USR100-43	4 x 3	96	140 x 566 x 318
		Total: 16 racks	-	(= 160)	-
ULUF P620MV	90/96/100	10 x USR100-42	4 x 2	80	140 x 566 x 212
		10 x USR100-43	4 x 3	120	140 x 566 x 318
		Total: 20 racks	-	(= 200)	-
ULUF P830MV	90/96/100	14 x USR100-42	4 x 2	112	140 x 566 x 212
		14 x USR100-43	4 x 3	168	140 x 566 x 318
		Total: 28 racks	-	(= 280)	-
ULUF 450-2M	90/96/100	15 x USR100-42	4 x 2	120	140 x 566 x 212
ULUF 750-2M	90/96/100	25 x USR100-42	4 x 2	200	140 x 566 x 212