

true cooling specialists



ULTRA LOW TEMPERATURE I CHEST FREEZER

Optimized internal storage capacity with stainless steel racks. Offers a lot of great controller features with the user-friendly controller in 4,3" display screen.

CONTROLLER FEATURES

Arctiko's ULTF series comes with the 4,3" G214 controller which combines security, monitoring & data logging functions.

This range allows data logging for up to 30 years, multiple alarms and much more customizability ensuring that the functionality of the chest models have the same features as the upright units in the same temperature range.

SUPERIOR COOLING PERFORMANCE

Reaching -80°C in only 190 minutes due to our our single compressor technology. The technology is specially designed to be the best in class in cooling performance.

With the ULTF freezer you are secured superior temperature stability and great temperature uniformity even after door openings. It also helps provide a stable cooling performance while lowering the energy consumption.

OPTIMIZED STORAGE CAPACITY

The compact space-saving cabinet design guarantees a reduced footprint because of the Arctiko's specifically developed vacuum insulated panels, optimizing the freezer to increase energy efficiency and cooling performance.

The vacuum insulated panels are of a high-efficient quality that secures the lowest possible energy consumption. VIP panels combined with Arctiko's state-of-the-art single compressor technology provide a safe and high-quality ULT freezer.




stored samples with

the stainless steel

front rack storage

system.

ULTE 220* V



ARCTIK CO ULTF 220 true cooling specialists


Pull-down times measures the time taken for temperature to rise, the faster the better. Hold over indication the freezers ability to keep the samples safe.

IMPROVED RECOVERY TIMES

.... Hold over

1110 1184 1258 1332 406

Pull-down



	MODEL	ULTF 220
Main specifications	Item code	DAI 0213
	Temperature range (°C)	-40 / -86
	Max. ambient temperature ("C)	20
	Capacity (L)	189
	Exterior dimensions WxDxH (mm)	920x691x890
	Interior dimensions WxDxH (mm)	760x440x630
	Weight (KG)	64
	Insulation (mm)	80
	Cooling technology	single system / static
	Controller model / probe type	G214 / PT1000
Power	Power supply (V)	230
	Frequency (Hz)	50/60
	Power consumption (kWh/24h)	5,9
	AMP (A)	2,7
	Power (Watt)	353
	Horse power (HP)	1/2
-	Alarm display as text - not codes	yes
Alarm functions	Visual / acoustic alarm	yes / yes
	Power failure alarm	yes
	Adjustable high / low temperature alarm	yes / yes
	Open door alarm	yes
	Probe failure alarm	yes
	Low battery alarm	yes
	Compressor failure alarm	yes
	Battery backup for alarms (approx. hrs.)	72
Data logging & external connection	Voltage free contact for remote alarm (e.g. GSM alarm module)	yes
	Electronic data logger	yes
	RS 485 port	yes
	USB port for software update & data download	yes
	Temperature chart recorder	optional
Display features	Display type	4,3" with touch buttons
	Number of optional reference probes	2
	Password protection for turning unit on / off	yes
	Temperature graph on display	yes
Other features	Automatic defrost	no
	Antifreeze function	no
	Vacuum release port	no
	Access port for external probe etc. (int. diameter mm)	12
	Light	по
	Lock with key	yes
	Castors (total pcs) / castors with brake (pcs)	4 / optional
		no
	Fixed feet (pcs) / levelling feet (pcs)	
	Sub lids / inner doors (pcs)	3
	Reversible door(s)	no
	Total capacity of 50 mm boxes 9x9 didviders (pcs)	117
	Total capacity of 75 mm boxes 9x9 didviders (pcs)	75
	Total capacity of 96 mm boxes 9x9 didviders (pcs)	54
Refrigeration & cabinet	Number of compressors	1
	Refrigerant	Nature R
	Interior cabinet material	painted steel
	Exterior cabinet material	painted steel
	Colour (exterior cabinet)	white
	Foam type	polyurethane
	VIP vacuum insulated panels	no
Shipping	Shipping weight (kg)	104
	Shipping dimensions WxDxH (cm)	102x76x106
	Shipping volume (m ³)	0,82

Arctiko A/S ■ Oddesundvej 39 · DK-6715 Esbjerg N Phone: +45 7020 0328 ■ info@arctiko.com ■ arctiko.com

10

-10

-20° -30°

-40

-50

-60

-70

-80 -90

4 148

- GSM alarm module
- Voltage converter
- Voltage stabilizer
- Racks & boxes
- Cryo gloves



CE

Actual performance of products may vary upon various factors. All rights to this temperature mapping belong to Arctiko A/S. The drawing may not be copied nor presented to any other person without prior written consent of Arctiko A/S.

CERTIFICATES