

true cooling specialists



## **ULTRA LOW TEMPERATURE I** CHEST FREEZER

Optimized internal storage capacity with stainless steel racks. Offers a lot of great controller features with the user-friendly controller in 4,3" display screen.

## **CONTROLLER FEATURES**

Arctiko's ULTF series comes with the 4,3" G214 controller which combines security, monitoring & data logging functions.

This range allows data logging for up to 30 years, multiple alarms and much more customizability ensuring that the functionality of the chest models have the same features as the upright units in the same temperature range.

## **SUPERIOR COOLING** PERFORMANCE

Reaching -80°C in only 190 minutes due to our our single compressor technology. The technology is specially designed to be the best in class in cooling performance.


With the ULTF freezer you are secured superior temperature stability and great temperature uniformity even after door openings. It also helps provide a stable cooling performance while lowering the energy consumption.

## OPTIMIZED STORAGE CAPACITY

The compact space-saving cabinet design guarantees a reduced footprint because of the Arctiko's specifically developed vacuum insulated panels, optimizing the freezer to increase energy efficiency and cooling performance.

The vacuum insulated panels are of a high-efficient quality that secures the lowest possible energy consumption. VIP panels combined with Arctiko's state-of-the-art single compressor technology provide a safe and high-quality ULT freezer.







NOISE REDUCTION TECHNOLOGY

User friendly G214 controller standard in ULTF range. Includes different alarms as; visual acoustic open door and more.




Lock on handle for optimal security of the valued samples. Comes as a standard feature.



Optimal storage capacity and overview of the stored samples with the stainless steel front rack storage system.





ARCTIK CO ULTF 420 true cooling specialists

148 222 296 370 592 592 666 666 666 814 814 818 888 888 888 8036

Pull-down times measures the time taken for temperature to rise, the faster the better. Hold over indication the freezers ability to keep the samples safe.

IMPROVED RECOVERY TIMES

Hold over

Pull-down

45

1110 1258 332 1406

- Racks & boxes
- Cryo gloves



CE

Arctiko A/S ■ Oddesundvej 39 · DK-6715 Esbjerg N Phone: +45 7020 0328 ■ info@arctiko.com ■ arctiko.com

	-75°	t																
()	-76°																	
	-77								ļ									
	-78°															4		
	-79°																	
Temperature ('C)	-80°	-	4	~	4	~	~		~	~	-	<u>/</u>	4	~		~	~	1
adma	-81°											ļ				-de-		{
4	-82°		·												tạt	pie	rur	1
	-83°																	{
	-84°		·			<u></u>	···	·		·	···		·					
	-85°	Ļ																
			28	58	88	118	148	178	208	238	268	298	328	358	388	418	448	
									Mi	nutes								

	MODEL	ULTF 420				
	Item code	DAI 0210				
	Temperature range (°C)	-40 / -86				
	Max. ambient temperature (°C)	20				
tions	Capacity (L)	368				
cifica	Exterior dimensions WxDxH (mm)	1562x691x890				
Main specifications	Interior dimensions WxDxH (mm)	1400x440x630				
	Weight (KG)	88				
	Insulation (mm)	80				
	Cooling technology	single system / static				
	Controller model / probe type	G214 / PT1000				
Power	Power supply (V)	230				
	Frequency (Hz)	50/60				
	Power consumption (kWh/24h)	8,5				
Ро	AMP (A)	2,6				
	Power (Watt)	490				
	Horse power (HP)	3/4				
	Alarm display as text - not codes	yes				
	Visual / acoustic alarm	yes / yes				
SL	Power failure alarm	yes				
lctio	Adjustable high / low temperature alarm	yes / yes				
Alarm functions	Open door alarm	yes				
	Probe failure alarm	yes				
	Low battery alarm	yes				
	Compressor failure alarm	yes				
	Battery backup for alarms (approx. hrs.)	72				
& ction	Voltage free contact for remote alarm (e.g. GSM alarm module)	yes				
ging	Electronic data logger	yes				
Data logging & external connection	RS 485 port	yes				
	USB port for software update & data download	yes				
ē	Temperature chart recorder	optional				
> s	Display type	4,3" with touch buttons				
Display features	Number of optional reference probes	2				
бD	Password protection for turning unit on / off	yes				
	Temperature graph on display	yes				
	Automatic defrost	no				
	Antifreeze function	no				
Other features	Vacuum release port Access port for external probe etc. (int. diameter mm)	no 12				
	Light	12 no				
	Lock with key					
	Castors (total pcs) / castors with brake (pcs)	yes 4 / optional				
	Fixed feet (pcs) / levelling feet (pcs)	no				
	Sub lids / inner doors (pcs)	4				
	Reversible door(s)	no				
	Total capacity of 50 mm boxes 9x9 didviders (pcs)	282				
	Total capacity of 75 mm boxes 9x9 didviders (pcs)	180				
	Total capacity of 96 mm boxes 9x9 didviders (pcs)	129				
Refrigeration & cabinet	Number of compressors	1				
	Refrigerant	Nature R				
	Interior cabinet material	painted steel				
	Exterior cabinet material	painted steel				
	Colour (exterior cabinet)	white				
	Foam type	polyurethane				
	VIP vacuum insulated panels	no				
Shipping	Shipping weight (kg)	132				
	Shipping dimensions WxDxH (cm)	166x76x106				
Shi	Shipping volume (m <sup>3</sup> )	1,34				

10' 0

-10 -20' -30' -40 -50' -60'

-70 -80' -90'

7

CERTIFICATES

Actual performance of products may vary upon various factors. All rights to this temperature mapping belong to Arctiko A/S. The drawing may not be copied nor presented to any other person without prior written consent of Arctiko A/S.